

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	1 of 32

Summary Of Safety and Clinical Performance Intended for Users

Reference

**EU MDR 2017/745, Article 32 - Summary of safety and clinical performance &
MDCG 2019-9 Rev.1 Summary of safety and clinical performance A guide for
manufacturers and notified bodies**

Product Details

Device Name	Variant Name
Endopyelotomy Stent- With/without Hydrophilic coated	Endopyelotomy Stent (External)
	Endoureterotomy Stent (Internal)

Product Classification

Class IIb, Rule 05 as per Annex VIII of MDR 2017/745

Manufacturer Details	Authorized Representative
1. Unit-I DEVON INNOVATIONS PRIVATE LIMITED No. 27A, Near State Bank of India, Electronic City Phase I, Hosur Main Road, Bangalore-560 100, India. Phone no: 080-28522354/28522367/28522368 2. Unit-II DEVON INNOVATIONS PRIVATE LIMITED Gupta complex, 1st floor, Khasra No: 519/370 Near EWS flats, sector-1, village Kamli Parwanoo 173220 Himachal Pradesh, India. Phone no: 01792232492	Amstermed BV Saturnusstraat 46-62, Unit 032,2132 HB Hoofddorp The Netherlands. Ph: +31 23 56 56 337 Email: info@amstermed.nl Website: https://www.amstermed.nl

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	2 of 32

Approvals	Name	Function	Signature	Date
Prepared By	Mrs. Janaki	QA/RA-Documentation in-charge		13.06.2025
Reviewed By	Mr. Srinivas	QA/RA-Manager		13.06.2025
Approved By	Mr. Ashwin Khemani	Director		13.06.2025

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
		Page No.:	3 of 32
	Summary of Safety and Clinical Performance		

Contents

Introduction.....	6
1. Device identification and general information	6
1.1 Device Trade Name(s)	6
1.2 Manufacturer's Name & Address	7
1.3 Manufacturer's Single Registration Number (SRN)	7
1.4 Basic UDI-DI	7
1.5 Medical Device Nomenclature	7
1.6 Class of device	7
1.7 Year of first certificate (CE) of the subject device	8
1.8 Authorized Representative	8
1.9 NB Details	8
1.10 Conformity Assessment Procedure	8
1.11 Link to SSCP in website	8
2. Intended use of the device.....	8
2.1 Intended Purpose	8
2.2 Indications & Target Populations	9
2.3 Contraindications	9
3. Device Description.....	9
3.1 Description of the Device	9
3.1.1 Device Models	10
3.1.2 Principle of Operation	10
3.2 Reference to previous generation(s) or variants.....	10
3.3 Accessories Details	11

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	4 of 32

3.4	Combination with other Medical Devices	11
4.	Risks and warnings	11
4.1	Residual risks and undesirable side effects	11
4.2	Warnings.....	12
4.3	Precautions	12
4.4	Other relevant aspects of safety	12
5	Summary of clinical evaluation and post-market clinical follow-up (PMCF)	13
5.1	Summary of clinical data related to similar device, if applicable	13
5.2	Summary Of Clinical Data from Conducted Investigations of the device before The CE-Marking 13	
5.3	Summary of clinical data from other sources	13
5.4	An overall summary of the clinical performance and safety.....	17
5.5.	Ongoing or planned post-market clinical follow-up.....	18
6.	Possible diagnostic or therapeutic alternative.....	20
7.	Suggested profile and training for users	20
8.	Reference to any harmonized standards and CS applied	20
8.1	Applicable Harmonized Standards	20
8.2	Other Applicable Standards.....	22
8.3	List of Guidelines	23
	Summary Of Safety and Clinical Performance	25
	Intended for Patients	25
	Introduction.....	25
1.	Device identification and general information	26
2.	Intended use of the device.....	26
2.1	Intended Purpose	26
2.2	Indications & Target Populations	26
2.3	Contraindications	27

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	5 of 32

3. Device Description	27
3.1 Device description	27
3.2 Materials that come in contact with patient	28
3.3 Information about medicinal substances in the device, if any	28
3.4 Description of how the device is achieving its intended mode of action	28
3.5 Description of accessories, if any	28
4. Risks and warnings	28
4.1 How potential risks have been controlled or managed	28
4.2 Residual Risks	29
4.3 Adverse events	29
4.4 Warnings	29
4.5 Precautions	30
4.6 Summary of any field safety corrective action, (FSCA including FSN) if applicable	30
5. Summary of clinical evaluation and post-market clinical follow-up	30
5.1 Clinical background of the device	30
5.2 The clinical evidence for the CE-marking	31
5.3 Safety	31
9. Possible diagnostic or therapeutic alternative	32
10. Suggested profile and training for users	32
11. Revision history	32

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	6 of 32

Introduction

The Regulation (EU) 2017/745 on medical devices requires that the manufacturer shall draw up a summary of safety and clinical performance (SSCP) for implantable devices and class III devices, other than custom-made or investigational devices. The SSCP shall be validated by a notified body (NB) and made available to the public via the European database on medical devices (Eudamed).

The SSCP is intended to provide public access to an updated summary of clinical data and other information about the safety and clinical performance of the medical device. The SSCP will be an important source of information for intended users – both healthcare professionals and relevant for patients. It is one of several means intended to fulfill the objectives of the Medical Device Regulation (MDR) to enhance transparency and provide adequate access to information.

The SSCP is not intended to:

- Give general advice on the diagnosis or treatment of particular medical conditions, nor
- Replace the instructions for use (IFU) as the main document that will be provided to ensure the safe use of a particular device, nor
- Replace the mandatory information on implant cards or in any other mandatory documents.

The main purpose of this document is to guide the presentation, content and validation of the SSCP. The word “shall” is used when there is a corresponding “shall” in the MDR, otherwise “should” or “recommended” etc. is used to indicate the interpretation of the MDR.

The following information is intended for users/healthcare professionals.

1. Device identification and general information

1.1 Device Trade Name(s)

Product Name	Endopyelotomy Stent- With/without Hydrophilic coated
Brand Name:	Devon
Variant Name:	Endopyelotomy Stent (External) and Endoureterotomy Stent (Internal)

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
		Page No.:	7 of 32
	Summary of Safety and Clinical Performance		

1.2 Manufacturer's Name & Address

Legal Manufacturer Name:	DEVON INNOVATIONS PRIVATE LIMITED
Registered Office & Manufacturing Unit-I Address:	No. 27A, Near State Bank of India, Electronic City Phase I, Hosur Main Road, Bangalore-560 100, India. Phone no: 080-28522354/28522367/28522368
Manufacturing Unit-II Address:	Gupta complex, 1st floor, Khasra No: 519/370, Near EWS flats, sector-1, village Kamli Parwanoo 173220 Himachal Pradesh, India. Phone no: 01792232492
Email:	srinivas@devoncath.com , nagendrakumar@devoncath.com
Website:	www.devoncath.com

1.3 Manufacturer's Single Registration Number (SRN)

Single Registration Number (SRN) for Unit-I Manufacturing site:	IN-MF-000010584
Single Registration Number (SRN) for Unit-II Manufacturing site:	IN-MF-000045808

1.4 Basic UDI-DI

8903410EPEX53

1.5 Medical Device Nomenclature

EMDN Code:	U020399
MDN Code:	MDN 1104-2
MDS and MDT code:	MDS 1005, MDT 2001, MDT 2002, MDT 2008, MDT 2011

1.6 Class of device

Class IIb, Rule 05, in accordance with Annex VIII of EU Medical Device Regulation 2017/745

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	8 of 32

1.7 Year of first certificate (CE) of the subject device

2012

1.8 Authorized Representative

Name:	Amstermed BV
Address:	Saturnusstraat 46-62, Unit 032, 2132 HB Hoofddorp The Netherlands.
Phone:	+31 23 56 56 337
Email:	info@amstermed.nl
Website:	https://www.amstermed.nl
SRN:	NL-AR-000001971

1.9 NB Details

Name:	DNV Product Assurance AS
Address:	Veritasveien 1, 1363 Høvik, Norway
Website:	www.dnv.com
Notified Body No.:	2460

1.10 Conformity Assessment Procedure

Conformity assessment procedure followed is Annexure IX.

1.11 Link to SSCP in website

The link for the Summary of Safety and Clinical Performance (SSCP) is provided below:

2. Intended use of the device

2.1 Intended Purpose

Used for temporary internal drainage from the Ureteropelvic junction to the bladder. In case of Hydrophilic coated, it is to improve the ease of insertion.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	9 of 32

2.2 Indications & Target Populations

Indication:

a. Endopyelotomy Stent (External):

- Extrinsic compression of ureter
- Ureteral incision
- Ureteropelvic junction incision
- Stricture dilatation

b. Endoureterotomy Stent (Internal):

Endoureterotomy Stent has been employed to relieve obstruction in a variety of benign, malignant, and post-traumatic conditions. These stents may be placed using endoscopic techniques. In case of Hydrophilic coated, it is to improve the ease of insertion.

Target patient population: Pediatrics and Adult

2.3 Contraindications

- Contraindicated surgical candidate
- Unexplained hematuria
- Unrepaired ureteral avulsion

3. Device Description


3.1 Description of the Device

A stent is a hollow tube that maintains patency until healing can take place or an obstruction is relieved.

Product Image:

Product Name	Variant Name	Packing	Image
Endopyelotomy Stent-With/without	Endopyelotomy Stent (External)	With Packing	

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
		Summary of Safety and Clinical Performance	Page No.:

Product Name	Variant Name	Packing	Image
Hydrophilic coated		Without Packing	
		With Packing	
	Endoureterotomy Stent (Internal)	Without Packing	

3.1.1 Device Models

#	Variant Name	Variant Description/Details
1.	Endopyelotomy Stent (External)	A stent is a hollow tube that maintains patency until healing can take place or an obstruction is relieved.
2.	Endoureterotomy Stent (Internal)	

3.1.2 Principle of Operation

A Endopyelotomy stent is a thin, flexible plastic tube that is inserted into the ureter to drain urine from the kidney into the bladder in the case of obstruction.

3.2 Reference to previous generation(s) or variants

Legacy Device Name:	Stents (Urology)- Endopyelotomy Stent- With/without Hydrophilic coated
Brand/Proprietary Name:	Devon
93/42/EEC (MDD) Cert. No.:	246182-2017-CE-IND-NA-PS, Rev.2.0
Notified Body Details:	DNV Product Assurance AS
Is any significant difference between Legacy Device &	There is no significant difference between legacy device and subject device with respect to raw materials used in production, device

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	11 of 32

Device Under Evaluation?

description, intended purpose, medical indications, target user, target patient population, side-effects, and contraindications.

3.3 Accessories Details

Not Applicable. The Endopyelotomy Stent- With/without Hydrophilic coated does not have any accessories supplied by the manufacturer.

3.4 Combination with other Medical Devices

Endopyelotomy Stent- With/without Hydrophilic coated is used along with Guidewire and positioning sleeve and guiding catheter. The stent can be used with any compatible guidewire and positioning sleeve and guiding catheter from other manufacturers.

4. Risks and warnings

4.1 Residual risks and undesirable side effects

a. Residual Risks

- Urinary Tract Infection
- Stent migration
- Contamination or Deterioration of product
- Toxic to environment
- Stent fragmentation
- Haematuria
- Tissue damage
- Delay in procedure, Inconvenience to the user

b. Adverse Events

- Migration
- Sepsis
- Encrustation

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	12 of 32

4.2 Warnings

All components of the Endopyelotomy Stent-With/ without Hydrophilic coated are for single use only.

- Duration of Use:
 - Periodic evaluation is advised. The Stent must not remain indwelling more than three months. These stents are not indented as permanent indwelling devices
 - Do not use device if there is any indication that the sterility of the device has been compromised.
- Adverse effects: Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance. Follow up procedures.
- Reuse: Reusing single-use stents can lead to urinary tract infections in patients.

4.3 Precautions

Carefully read all instructions for use and product labeling. The device shall only be applied for its intended use, and in accordance with these instructions. Observe all cautions and warnings throughout these instructions. Failure to do so may result in complications.

All Health care professionals is responsible for using the appropriate technique and deciding on the indication for use of this device based on own experience, training and medical judgment. The doctor must be trained in the proper use of the device.

4.4 Other relevant aspects of safety

There were no identified and/or received reportable events that led to death, a serious deterioration in the state of health of the patient, user, or other person for Endopyelotomy Stent-With/ without Hydrophilic coated. Hence FSCA or FSN is not applicable.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	13 of 32

5 Summary of clinical evaluation and post-market clinical follow-up (PMCF)

5.1 Summary of clinical data related to similar device, if applicable

The Endopyelotomy Stent-With/ without Hydrophilic coated belongs to the “Urology” group. In the present market there are many similar devices and/or benchmark devices available with same intended purpose and are having the same generally acknowledged state-of-the-art.

These similar devices fall under “Well-Established Technology”. Data from similar devices is considered for the conformation of conformity to the Endopyelotomy Stent-With/ without Hydrophilic coated relevant general safety and performance requirements. The similar device data is used to demonstrate ubiquity of design, lack of novelty, known safety and performance profile of a generic group of devices, etc.

The below mentioned similar devices data is used to evaluate the Endopyelotomy Stent-With/ without Hydrophilic coated relevant general safety and performance requirements as part of literature review. These similar devices contain the same Raw materials and same intended purpose, but due to insufficient information availability the full assessment of equivalence is not possible. Therefore, these similar devices are used for the same clinical intended purposes as the Endopyelotomy Stent-With/ without Hydrophilic coated and are considered to be similar but non-equivalent devices.

S.NO	Product Name	Similar device	Manufacturer Name
1.	Endopyelotomy Stent (External)	Endopyelotomy Stent	Boston Scientific Corporation
2.	Endoureterotomy Stent (Internal)	Endoureterotomy Stent	Cook Incorporated

The similar device’s SSCP would be available in EUDAMED.

5.2 Summary Of Clinical Data from Conducted Investigations of the device before The CE-Marking

Not Applicable

5.3 Summary of clinical data from other sources

The below mentioned are literatures selected for detailed review for:

- Evaluation of state of the art
- Evaluation of clinical data from similar devices

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	14 of 32

A. Endopyelotomy Stent (External)

#	ID#	Source Link	Literature Title
1.	AL1	https://pubmed.ncbi.nlm.nih.gov/15245924/	Endopyeloplasty Versus Endopyelotomy Versus Laparoscopic Pyeloplasty for Primary Ureteropelvic Junction Obstruction
2.	AL3	https://pubmed.ncbi.nlm.nih.gov/11490243/	Influence Of Stent Size on The Success of Antegrade Endopyelotomy for Primary Ureteropelvic Junction Obstruction: Results Of 2 Consecutive Series
3.	AL5	https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3516940/pdf/WIITM-6-17369.pdf	Ureteroscopic Holmium Laser Endopyelotomy for Ureteropelvic Junction Stenosis After Pyeloplasty
4.	AL7	https://link.springer.com/article/10.1007/s11934-007-0062-0	Antegrade Percutaneous Endopyelotomy
5.	AL9	https://www.sciencedirect.com/science/article/pii/S0090429503002310	Impact Of Hydronephrosis and Renal Function on Treatment Outcome: Antegrade Versus Retrograde Endopyelotomy
6.	AL10	https://www.sciencedirect.com/science/article/pii/S0090429505009714	Laparoscopic Pyeloplasty Versus Antegrade Endopyelotomy: Comparison In 100 Patients and A New Algorithm for The Minimally Invasive Treatment of Ureteropelvic Junction Obstruction
7.	AL11	https://pubmed.ncbi.nlm.nih.gov/11697405/	Endopyelotomy in Poorly Functioning Kidney: Is It Worthwhile?
8.	AL12	https://www.sciencedirect.com/science/article/pii/S0377123703801454	Endopyelotomy - a Minimally Invasive Surgical Option for Pelvi-ureteric Junction Obstruction: a Study Of 34 Cases
9.	AL13	https://pubmed.ncbi.nlm.nih.gov/12576806/	Is A 2-Week Duration Sufficient for Stenting In Endopyelotomy?
10.	AL14	https://pubmed.ncbi.nlm.nih.gov/15801361/	Endopyelotomy in Childhood: Our Experience with 37 Patients
11.	AL15	https://www.sciencedirect.com/science/article/pii/S0090429500007585	Simplified Approach to Percutaneous Endopyelotomy

	Devon Innovations Private Limited		Document No.:	TD/DIP/SSCP/03
			Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated		Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance		Page No.:	15 of 32

#	ID#	Source Link	Literature Title
12.	AL16	https://www.goldjournal.net/article/S0090-4295(04)00321-8/abstract	Update on ureteral stents
13.	AL17	https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5734136/	Massive Dilation of the Ureter: An Endoscopic Management of Persistent Urinary Leak After Partial Nephrectomy
14.	AL18	https://pubmed.ncbi.nlm.nih.gov/25246158/	Laser endoureterotomy and endopyelotomy: an update
15.	AL19	https://www.researchgate.net/publication/381365518_Contemporary_status_of_diagnostic_endoluminal_ultrasound_and_optical_coherence_tomography_in_the_ureter	Contemporary status of diagnostic endoluminal ultrasound and optical coherence tomography in the ureter
16.	AL20	https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4198764/	Long-Term Follow-up Results of Laparoscopic Pyeloplasty

B. Endoureterotomy Stent (Internal)

#	ID#	Source Link	Literature Title
1.	BL1	https://pubmed.ncbi.nlm.nih.gov/26102617/	Single Versus Double Ureteral Stent Placement After Laser Endoureterotomy for the Management of Benign Ureteral Strictures: A Randomized Clinical Trial
2.	BL3	https://pubmed.ncbi.nlm.nih.gov/15720342/	Ureteroscopic endoureterotomy
3.	BL5	https://pubmed.ncbi.nlm.nih.gov/24294293/	Current status of minimally invasive endoscopic management of ureteric strictures
4.	BL7	https://pubmed.ncbi.nlm.nih.gov/12376232/	Endoscopic Treatment of Benign Ureteral Strictures
5.	BL9	https://pubmed.ncbi.nlm.nih.gov/10795616/	Endoureterotomy for Congenital Primary Obstructive Megaureter: Preliminary Report.
6.	BL11	https://pubmed.ncbi.nlm.nih.gov/17705763/	Endoureterotomy for Treatment of Primary Obstructive Megaureter in Children

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	16 of 32

#	ID#	Source Link	Literature Title
7.	BL12	https://pubmed.ncbi.nlm.nih.gov/23601441/	Endoureterotomy as the Initial Management of Concurrent Ureteropelvic and Ureterovesical Junction Obstruction After Failed Conservative Therapy
8.	BL13	https://pubmed.ncbi.nlm.nih.gov/29500521/	Endoscopic Management of Ureteral Strictures: An Update
9.	BL14	https://pubmed.ncbi.nlm.nih.gov/25246158/	Laser endoureterotomy and endopyelotomy: an update
10.	BL15	https://pubmed.ncbi.nlm.nih.gov/15245923/	Update On Ureteral Stents
11.	BL16	https://pubmed.ncbi.nlm.nih.gov/PMC7292192/	Dual ureteral stent placement after redo laser endoureterotomy to manage persistent ureteral stricture
12.	BL17	https://pubmed.ncbi.nlm.nih.gov/36246791/	The “Cut-to-the-Light” Technique Laser Endoureterotomy for Complete Ureteral Obstruction Resurfaces! A New Application of an Old Technique
13.	BL18	https://pubmed.ncbi.nlm.nih.gov/19419288/	Ureteral stent: past, present and future
14.	BL19	https://pubmed.ncbi.nlm.nih.gov/31549458/	Ureteric stents: Overview of current clinical applications and economic implications
15.	BL20	https://www.frontiersin.org/articles/10.3389/fruro.2023.1150795/full	Endourological treatment of upper tract urinary disease in children

Examples of clinical data registries are:

1. Clinicaltrials.gov: ClinicalTrials.gov is a database of privately and publicly funded clinical studies conducted around the world. ClinicalTrials.gov is a resource provided by the U.S. National Library of Medicine.
2. Clinicaltrialsregister.eu: The European Union Clinical Trials Register allows you to search for protocol and results information on interventional clinical trials that are conducted in the European Union (EU) and the European Economic Area (EEA); clinical trials conducted outside the EU / EEA that are linked to European pediatric-medicine development.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	17 of 32

3. Ctri.gov: The Clinical Trials Registry- India (CTRI), hosted at the ICMR's National Institute of Medical Statistics (<http://icmr-nims.nic.in>), is a free and online public record system for registration of clinical trials being conducted in India. Any researcher who plans to conduct a trial involving human participants, of any intervention such as drugs, surgical procedures, preventive measures, lifestyle modifications, devices, educational or behavioral treatment, rehabilitation strategies as well as trials are expected to register the trial in the CTRI before enrollment of the first participant. In the CTRI, details of Indian investigators, trial sites, Indian target sample size and date of enrollment are captured.

We have not found any studies related to the product.

The advantage of referring the data from a third-party registry is usually centered on patients or cases, which provides a sufficient level of traceability and detail for the study of our product. Using registries can be less expensive than initiating new and private device registries. The data in these registries might be relevant but may not contain the exact attributes needed to document/demonstrate clinical performance.

Hence, we cannot completely rely on the above registries to collect enough data for our devices. We have referred these registries only to gather information on similar devices and the study methodology.

5.4 An overall summary of the clinical performance and safety

Residual Risks	Medical Benefits
1. Urinary Tract Infection 2. Stent migration 3. Contamination or Deterioration of product 4. Toxic to environment 5. Stent fragmentation 6. Haematuria 7. Tissue damage 8. Delay in procedure, Inconvenience to the user	Subject Device 1. Relief of urinary obstruction 2. Continuous urine drainage from kidney 3. Prevention of Ureteral Stricture Formation Similar Device 1. Shorter operative time 2. Decreased morbidity 3. Safe and effective procedure for UPJ obstruction 4. Minimally invasive treatment 5. Long-term success rate 6. Symptomatic improvement 7. Improved drainage 8. Reduced postoperative pain 9. Faster convalescence Improvement of GFR (Glomerular filtration rate)

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	18 of 32

All the residual risks were reviewed and analysed and also identified the medical benefits of the intended use outweigh the overall residual risk. Endopyelotomy Stent- With/without Hydrophilic coated complies with all the Safety and Performance requirements with respect to the intended use of the device from the Clinical Evaluation study. The clinical evaluation is complete and conforms to the general safety and performance requirements. The Clinical evidence is demonstrated with the relevant General Safety & Performance Requirements as per Annex I of MDR. Risk Mitigation has been established as per the guidelines of ISO 14971.

5.5. Ongoing or planned post-market clinical follow-up

Overall, 50 subjects have participated in this PMCF study of the product Endopyelotomy Stent- With/without Hydrophilic coated. As per the PMCF study, the Endopyelotomy Stent- With/without Hydrophilic coated have met the primary and secondary objectives.

Overall study results:

Parameter	Study Results																					
Subjects	A total of 50 Subjects were enrolled for the PMCF study. As per the study, 43 (86%) subjects were male and 07 (14%) subjects were female.																					
	Subjects Group summary – 38 (76%) subjects were adults and 12 (24%) subjects were Pediatrics.																					
	Age Group:																					
	<table><tr><th>Age Group</th><th>No: of subjects</th><th>Percentage (%)</th></tr><tr><td>11-18 Years</td><td>12</td><td>24%</td></tr><tr><td>19-30 years</td><td>0</td><td>0</td></tr><tr><td>31-40 Years</td><td>09</td><td>18%</td></tr><tr><td>41-50 Years</td><td>24</td><td>48%</td></tr><tr><td>51-60 Years</td><td>04</td><td>8%</td></tr><tr><td>61-70 Years</td><td>01</td><td>2%</td></tr></table>	Age Group	No: of subjects	Percentage (%)	11-18 Years	12	24%	19-30 years	0	0	31-40 Years	09	18%	41-50 Years	24	48%	51-60 Years	04	8%	61-70 Years	01	2%
	Age Group	No: of subjects	Percentage (%)																			
	11-18 Years	12	24%																			
	19-30 years	0	0																			
	31-40 Years	09	18%																			
41-50 Years	24	48%																				
51-60 Years	04	8%																				
61-70 Years	01	2%																				
Target Users	Urologist																					
Study Site	National Kidney hospital, Jalandhar Yashoda hospital,Pune Kshotri Nursing Home, Pune																					

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	19 of 32

Clinical Indication	Clinical Indication of the subjects for using Endopyelotomy Stent/Set - With/without Hydrophilic Coated: - Extrinsic compression of ureter - Ureteral incision - Ureteropelvic junction incision - Stricture dilatation Clinical Indication of the subject for using Endoureterotomy Stent - To receive obstruction in a variety of benign, malignant, and post-traumatic conditions
Clinical Safety	All the subjects considered for this study benefitted out of using Endopyelotomy Stent- With/without Hydrophilic coated. This clinically proves the safety of using our product on subjects with better efficacy.
Clinical Performance	The overall rating was 4 which is “Good” as per the definition, hence the product performance was clinically proven to be “Good”.
Adverse events and Risks	No users reported any of the adverse events during the study. No reports on the risks imposed by our Endopyelotomy Stent- With/without Hydrophilic coated.
Follow Up Summary	The stent was removed during a follow up procedure. The subject had no issues with stent removal. The Endopyelotomy Stent- With/without Hydrophilic coated effectively relieved the obstructive symptoms, and the patient experienced significant symptom improvement after the procedure. No complications were observed during follow-up.
Product Experience	The overall rating is 8, it is “Good” as per the definition. Hence this proves that our users were satisfied with the product and its purpose.

The Endopyelotomy Stent- With/without Hydrophilic coated from Devon Innovations Private Limited has reached all the safety and performance requirements with respect to the intended use of the device from the Post Market Follow up Clinical study. The Performance and Safety of the device as we claimed have been established in the Technical File and Instruction for Use (IFU). There were no new risks identified from the PMCF study for the product hence there is no addition to the residual risks which we have already identified in the Risk Management Report and that is been mitigated and are acceptable when weighed against the benefits to the patient.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	20 of 32

6. Possible diagnostic or therapeutic alternative

a. Endopyelotomy Stent (External):

- Open Pyeloplasty
- Laproscopic Pyeloplasty
- Retrograde Endopyelotomy
- Percutaneous Nephrostomy

b. Endoureterotomy Stent (Internal):

- Open or laparoscopic surgical repair
- Balloon dilation alone
- Endoureterotomy without stent placement

7. Suggested profile and training for users

Target Users: Urologist

Urologists diagnose and treat diseases of the urinary tract in both men and women. They also diagnose and treat anything involving the reproductive tract in men.

The target users are aware of basic operations of Endopyelotomy Stent- With/without Hydrophilic coated. There is no special user training is required. However, the device related directions for use information are provided in the Instruction for Use.

8. Reference to any harmonized standards and CS applied

8.1 Applicable Harmonized Standards

#	Standard ID	Current Issue	Title
Quality Management System Requirements			
1.	EN ISO 13485	2016/AC:2018/A11:2021	Medical devices - Quality management systems - Requirements for regulatory purposes (ISO 13485:2016)
Risk Management Requirements			
2.	EN ISO 14971	2019/A11:2021	Medical devices - Application of risk management to medical devices (ISO

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
		Summary of Safety and Clinical Performance	Page No.:

#	Standard ID	Current Issue	Title
			14791:2019)
Biological Risk Evaluation Requirements			
3.	EN ISO 10993-10	2023	Biological evaluation of medical devices - Part 10: Tests for skin sensitization (ISO 10993-10:2021)
4.	EN ISO 10993-23	2021	Biological evaluation of medical devices - Part 23: Tests for irritation (ISO 10993-23:2021)
Labels & Symbols Requirements			
5.	EN ISO 15223-1	2021	Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
Packaging Requirements			
6.	EN ISO 11607-1	2020/A1:2023	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems - (ISO 11607-1:2019)
7.	EN ISO 11607-2	2020/A1:2023	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes (ISO 11607-2:2019)
EO Sterilization Requirements			
8.	EN ISO 11135	2014/A1:2019	Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices (ISO11135:2014)
Sterility Test Requirements			
9.	EN ISO 11737-1	2018/A1:2021	Sterilization of health care products - Microbiological methods - Part 1: Determination of a population of microorganisms on products - (ISO 11737-1:2018)
10.	EN ISO 11737-2	2020	Sterilization of health care products - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process (ISO 11737-2:2019)

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	22 of 32

8.2 Other Applicable Standards

#	Standard ID	Current Issue	Title
Risk Management Requirements			
1.	ISO/TR 24971	2020	Medical devices – Guidance on the application of ISO 14971 (ISO/TR 24971:2020)
Usability			
2.	EN 62366-1	2015/A1:2020	Medical devices - Part 1: Application of usability engineering to medical devices
3.	IEC 62366-1	2015/A1:2020	Medical devices — Part 1: Application of usability engineering to medical devices Amendment 1
Biological Risk Evaluation Requirements			
4.	EN ISO 10993-1	2020	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process (ISO 10993-1:2018, including corrected version 2018-10)
5.	EN ISO 10993-3	2014	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity (ISO 10993-3:2014)
6.	EN ISO 10993-5	2009/A11:2025	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
7.	EN ISO 10993-6	2016	Biological evaluation of medical devices - Part 6: Tests for local effects after implantation (ISO 10993-6:2016)
8.	EN ISO 10993-7	2008/A1:2022	Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals - Amendment 1: Applicability of allowable limits for neonates and infants (ISO 10993-7:2008/Amd 1:2019)
9.	EN ISO 10993-11	2018	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity (ISO 10993-11:2017)
Instructions For Use Requirements			
10.	EN ISO 20417	2021	Medical devices - Information to be supplied by the manufacturer (ISO 20417:2021, Corrected version 2021-12)
Medical Device “Sterile” Requirements			
11.	EN 556-1	2024	Sterilization of medical devices - Requirements

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	23 of 32

#	Standard ID	Current Issue	Title
			for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
Transport Requirements			
12.	ISTA 2A	2011	Partial Stimulation Performance Test Procedure Packaged Products weighing 150 lbs (68 kg) or less
Cleanroom Requirements			
13.	ISO 14644-1	2015	Clean Rooms and Associated Controlled Environments – Part 1: Classification of Air Cleanliness by Particle Concentration
14.	ISO 14644-2	2015	Clean Rooms and Associated Controlled Environments – Part 2: Monitoring to Provide Evidence of Clean Room Performance Related to Air Cleanliness by Particle Concentration
Post Market Surveillance Requirements			
15.	ISO/ TR 20416	2020	Medical devices – Post market surveillance for manufacturers
Stability Requirements			
16.	ASTM F 1980-21	2021	Standard Guide for Accelerated Aging of Sterile Barrier Systems and Medical Devices
Labels & Symbols Requirements			
17.	ISO 15223-1	2021/Amd 1:2025	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements Amendment 1: Addition of defined term for authorized representative and modified EC REP symbol to not be country or region specific

8.3 List of Guidelines

#	Guideline	Current Issue	Title
1.	MEDDEV 2.7/1 Rev. 4	June 2016	Clinical Evaluation: A guide for manufacturers and notified bodies under directives 93/42/EEC and 90/385/EEC
2.	MEDDEV 2.5/5 Rev. 3	February 1998	Translation Procedure - Guidelines relating to the application of: The council directive 90/385/EEC on active implantable medical



Devon Innovations Private Limited

Endopyelotomy Stent- With/without Hydrophilic coated

Summary of Safety and Clinical Performance

Document No.:	TD/DIP/SSCP/03
Revision No.:	00
Effective Date:	13.06.2025
Page No.:	24 of 32

#	Guideline	Current Issue	Title
			devices The council directive 93/42/EEC on medical devices
3.	MEDDEV 2.12-1 Rev. 8	January 2013	Guidelines on a Medical Devices Vigilance System
4.	NB-MED 2.12/Rec. 1	February 2000	Post-Marketing Surveillance (PMS) post market/production
5.	MDCG 2021-24	October 2021	Guidance on classification of medical devices
6.	MDCG 2020-6	April 2020	Regulation (EU) 2017/745: Clinical evidence needed for medical devices previously CE marked under Directives 93/42/EEC or 90/385/EEC - A guide for manufacturers and notified bodies
7.	MDCG 2018-1 Rev.4	April 2021	Guidance on BASIC UDI-DI and changes to UDI-DI
8.	MDCG 2020-7	April 2020	Post-market clinical follow-up (PMCF) Plan Template - A guide for manufacturers and notified bodies
9.	MDCG 2020-8	April 2020	Post-market clinical follow-up (PMCF) Evaluation Report Template - A guide for manufacturers and notified bodies
10.	MDCG 2022-21	December 2022	Guidance on Periodic Safety Update Report (PSUR) according to Regulation (EU) 2017/745 (MDR)
11.	MDCG 2019-9 Rev.1	March 2022	Summary of safety and clinical performance A guide for manufacturers and notified bodies
12.	MDCG 2020-3 Rev.1 6	September 2023	Guidance on significant changes regarding the transitional provision under Article 120 of the MDR - May 2023
13.	MDCG 2024-2, Rev.1	January 2025	Procedures for the updates of the European Medical Device Nomenclature

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	25 of 32

Summary Of Safety and Clinical Performance

Intended for Patients

Introduction

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of its safety and clinical performance is prepared for users/healthcare professionals.

The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your healthcare professional in case you have questions about your medical condition or about the use of the device in your situation. This SSCP is not intended to replace an Implant card or the Instructions for Use to provide information on the safe use of the device.

The following information is intended for patient.

READABILITY

The SSCP has one part for intended users/healthcare professionals, and a second part for patients. It provides clear information at an appropriate depth to reflect the healthcare professionals' and the patients' different levels of knowledge. The readability of the part of the SSCP intended for patients is assessed by a test given to lay persons.

The readability was assessed by providing the document to laypersons with different literacy background to read. All the laypersons who have read the document were able to understand the terms and details of the device properly. Also, they agreed that the document language was legible and understandable. The manufacturer used this method to confirm that the SSCP was written in a way that is clear to the patient.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	26 of 32

1. Device identification and general information

Product Name	Endopyelotomy Stent- With/without Hydrophilic coated
Brand Name:	Devon
Models:	Endopyelotomy Stent (External) and Endoureterotomy Stent (Internal)
Manufacturer's Name & Address	DEVON INNOVATIONS PRIVATE LIMITED Registered Office & Manufacturing Unit-I Address: No. 27A, Near State Bank of India, Electronic City Phase I, Hosur Main Road, Bangalore-560 100, India. Phone no: 080-28522354/28522367/28522368 Manufacturing Unit-II Address: Gupta complex, 1st floor, Khasra No: 519/370, Near EWS flats, sector-1, village Kamli Parwanoo 173220 Himachal Pradesh, India. Phone no: 01792232492 Email: srinivas@devoncath.com nagendrakumar@devoncath.com Website: www.devoncath.com
Basic UDI-DI	8903410EPEX53
Year when the device was first CE-marked	2012

2. Intended use of the device

2.1 Intended Purpose

Used for temporary internal drainage from the Ureteropelvic junction to the bladder. In case of Hydrophilic coated, it is to improve the ease of insertion.

2.2 Indications & Target Populations

Indication:

- a. **Endopyelotomy Stent (External):**
 - Extrinsic compression of ureter
 - Ureteral incision

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	27 of 32

- Ureteropelvic junction incision
- Stricture dilatation

b. Endoureterotomy Stent (Internal):

Endoureterotomy Stent has been employed to relieve obstruction in a variety of benign, malignant, and post-traumatic conditions. These stents maybe placed using endoscopic, techniques. In case of Hydrophilic coated, it is to improve the ease of insertion.

Target patient population: Pediatrics and Adult

2.3 Contraindications


- Contraindicated surgical candidate
- Unexplained hematuria
- Unrepaired ureteral avulsion

3. Device Description

3.1 Device description

A stent is a hollow tube that maintains patency until healing can take place or an obstruction is relieved.

Product Image:

Product Name	Variant Name	Packing	Image
Endopyelotomy Stent-With/without Hydrophilic coated	Endopyelotomy Stent (External)	With Packing	
		Without Packing	
	Endoureterotomy Stent (Internal)	With Packing	

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
		Summary of Safety and Clinical Performance	Page No.:

Product Name	Variant Name	Packing	Image
		Without Packing	

3.2 Materials that come in contact with patient

Material that comes in contact with patient is stent and Guidewire.

3.3 Information about medicinal substances in the device, if any

Endopyelotomy Stent- With/without Hydrophilic coated does not incorporate medicinal substances. Hence this declaration is not applicable

3.4 Description of how the device is achieving its intended mode of action

Endopyelotomy is a surgical procedure aimed at treating pelvic-ureteric junction (PUJ) obstruction, which is a blockage where the renal pelvis meets the ureter. After an endopyelotomy procedure, endopyelotomy stents are often used to maintain the patency of the surgical site, promote healing, and facilitate urine flow from the kidney to the bladder.

3.5 Description of accessories, if any

Not Applicable. The Endopyelotomy Stent- With/without Hydrophilic coated does not have any accessories supplied by the manufacturer

4. Risks and warnings

Contact your healthcare professional if you believe that you are experiencing side-effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

4.1 How potential risks have been controlled or managed

We have established and maintaining a strong quality management system and continuous process monitoring controls as per EN ISO 13485:2016/A11:2021 requirements. As per the risk management

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	29 of 32

plan, all the necessary and possible risk control measures are implemented to reduce the risk to practicable acceptable level. We have implemented all necessary the risk control measures for each hazard for reducing the risks to an acceptable level.

The Endopyelotomy Stent- With/without Hydrophilic coated risk management is completed by risk estimation, risk analysis, risk evaluation, risk control, overall residual risk evaluation, production, and postproduction information. The Endopyelotomy Stent- With/without Hydrophilic coated is reached all the safety and performance claims when used as per the defined indented purpose.

The Endopyelotomy Stent- With/without Hydrophilic coated medical benefits and residual risks are compared based on medical condition, literature data certainty, similar device data (benefits & residual risks) from the market, acknowledged state of the art.

4.2 Residual Risks

- Urinary Tract Infection
- Stent migration
- Contamination or Deterioration of product
- Toxic to environment
- Stent fragmentation
- Haematuria
- Tissue damage
- Delay in procedure, Inconvenience to the user

4.3 Adverse events

- Migration
- Sepsis
- Encrustation

4.4 Warnings

All components of the Endopyelotomy Stent- With/without Hydrophilic coated are for single use only.

- Duration of Use:
 - Periodic evaluation is advised. The Stent must not remain indwelling more than three months. These stents are not indented as permanent indwelling devices
 - Do not use device if there is any indication that the sterility of the device has been compromised.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	30 of 32

- Adverse effects: Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance. Follow up procedures.
- Reuse: Reusing single-use stents can lead to urinary tract infections in patients.

4.5 Precautions

Carefully read all instructions for use and product labeling. The device shall only be applied for its intended use, and in accordance with these instructions. Observe all cautions and warnings throughout these instructions. Failure to do so may result in complications.

All Health care professionals is responsible for using the appropriate technique and deciding on the indication for use of this device based on own experience, training and medical judgment. The doctor must be trained in the proper use of the device.

4.6 Summary of any field safety corrective action, (FSCA including FSN) if applicable

There were no identified and/or received reportable events that led to death, a serious deterioration in the state of health of the patient or user, for Endopyelotomy Stent- With/without Hydrophilic coated. Hence FSCA or FSN is not applicable.

5. Summary of clinical evaluation and post-market clinical follow-up

5.1 Clinical background of the device

The Endopyelotomy Stent- With/without Hydrophilic coated is a designed and developed as per the latest and/or current technical, international and regulatory. The Endopyelotomy Stent- With/without Hydrophilic coated performance complies all necessary requirements required for its intended purpose. The Endopyelotomy Stent- With/without Hydrophilic coated related all known foreseeable hazards are identified and associated risk are reduced as far as possible by implementing all necessary risk control measures as per the requirements of EN ISO 14971:2019/A11:2021. The Endopyelotomy Stent- With/without Hydrophilic coated is manufactured as per the defined standard operating procedures in controlled environments by training personnel and complies all necessary requirements of EN ISO 13485:2016/A11:2021.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
		Page No.:	31 of 32
	Summary of Safety and Clinical Performance		

5.2 The clinical evidence for the CE-marking

Endopyelotomy Stent- With/without Hydrophilic coated is having a CE certified medical device under EU MDD 93/42/EEC. Hence, the Endopyelotomy Stent- With/without Hydrophilic coated comply the definition as a legacy device as per the MDCG 2020-6:2020 – Regulation (EU) 2017/745: Clinical evidence needed for medical devices previously CE marked under Directives 93/42/EEC or 90/385/EEC.

Device Name:	Stents (Urology)- Endopyelotomy Stent- With/without Hydrophilic coated
Brand/Proprietary Name:	Devon
93/42/EEC (MDD) Cert. No.:	246182-2017-CE-IND-NA-PS, Rev.2.0
Notified Body Details:	DNV Product Assurance AS

The Endopyelotomy Stent- With/without Hydrophilic coated belongs to the “Urology” group. In the present market there are many similar devices and/or benchmark devices available with same intended purpose and are having the same generally acknowledged state-of-the-art.

S.NO	Product Name	Similar device	Manufacturer Name
1.	Endopyelotomy Stent (External)	Endopyelotomy Stent	Boston Scientific Corporation
2.	Endoureterotomy Stent (Internal)	Endoureterotomy Stent	Cook Incorporated

5.3 Safety

We have reviewed and analysed all the residual risks and also identified the medical benefits of the intended use outweigh the overall residual risk. All the residual risks are acceptable by providing appropriate information to the end user’s awareness in the form of “Label’ and ‘Instructions for Use”.

Residual Risks	Medical Benefits
1. Urinary Tract Infection 2. Stent migration 3. Contamination or Deterioration of product 4. Toxic to environment 5. Stent fragmentation 6. Haematuria 7. Tissue damage 8. Delay in procedure,	Subject Device 1. Relief of urinary obstruction 2. Continuous urine drainage from kidney 3. Prevention of Ureteral Stricture Formation Similar Device 1. Shorter operative time 2. Decreased morbidity 3. Safe and effective procedure for UPJ obstruction 4. Minimally invasive treatment

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/03
		Revision No.:	00
	Endopyelotomy Stent- With/without Hydrophilic coated	Effective Date:	13.06.2025
	Summary of Safety and Clinical Performance	Page No.:	32 of 32

Residual Risks	Medical Benefits
Inconvenience to the user	5. Long-term success rate 6. Symptomatic improvement 7. Improved drainage 8. Reduced postoperative pain 9. Faster convalescence 10. Improvement of GFR (Glomerular filtration rate)

9. Possible diagnostic or therapeutic alternative

When considering alternative treatments, it is recommended to contact your healthcare professional who can take into account your individual situation.

10. Suggested profile and training for users

Target Users: Urologist

Urologists diagnose and treat diseases of the urinary tract in both men and women. They also diagnose and treat anything involving the reproductive tract in men.

The target users are aware of basic operations of Endopyelotomy Stent- With/without Hydrophilic coated. There is no special user training is required. However, the device related directions for use information are provided in the Instruction for Use.

11. Revision history

SSCP Rev. No.	Date Issued	Change description	Rev. Validated by the NB
00	13.06.2025	Initial Release	☐ Yes Validation language: English ☒ No