

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	1 of 31

Summary Of Safety and Clinical Performance Intended for Users

Reference

EU MDR 2017/745, Article 32 - Summary of safety and clinical performance & MDCG 2019-9 Rev.1 Summary of safety and clinical performance A guide for manufacturers and notified bodies

Product Details

Device Name	Variant Name
Urethral Stents- With/without Hydrophilic coated	-----

Product Classification

Class IIb, Rule 05 as per Annex VIII of MDR 2017/745

Manufacturer Details	Authorized Representative
1. Unit-I DEVON INNOVATIONS PRIVATE LIMITED No. 27A, Near State Bank of India, Electronic City Phase I, Hosur Main Road, Bangalore-560 100, India. Phone no: 080-28522354/28522367/28522368 2. Unit-II DEVON INNOVATIONS PRIVATE LIMITED Gupta complex, 1st floor, Khasra No: 519/370 Near EWS flats, sector-1, village Kamli Parwanoo 173220 Himachal Pradesh, India. Phone no: 01792232492	Amstermed BV Saturnusstraat 46-62, Unit 032,2132 HB Hoofddorp The Netherlands. Ph: +31 23 56 56 337 Email: info@amstermed.nl Website: https://www.amstermed.nl



**Devon Innovations Private
Limited**

**Urethral Stents-With/without
Hydrophilic coated**


**Summary of Safety and Clinical
Performance**

Document No.: TD/DIP/SSCP/02

Revision No.: 00

Effective Date: 23.06.2025

Page No.: **2 of 31**

Approvals	Name	Function	Signature	Date
Prepared By	Mrs. Janaki	QA/RA- Documentation in-charge		23.06.2025
Reviewed By	Mr. Srinivas	QA/RA-Manager		23.06.2025
Approved By	Mr. Ashwin Khemani	Director		23.06.2025

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	3 of 31

Contents

Introduction	6
1. Device identification and general information	6
1.1 Device Trade Name(s)	6
1.2 Manufacturer's Name & Address	7
1.3 Manufacturer's Single Registration Number (SRN)	7
1.4 Basic UDI-DI	7
1.5 Medical Device Nomenclature	7
1.6 Class of device	7
1.7 Year of first certificate (CE) of the subject device	8
1.8 Authorized Representative	8
1.9 NB Details	8
1.10 Conformity Assessment Procedure	8
1.11 Link to SSCP in website	8
2. Intended use of the device	8
2.1 Intended Purpose	8
2.2 Indications & Target Populations	9
2.3 Contraindications	9
3. Device Description	9
3.1 Description of the Device	9
3.1.1 Device Models	10
3.1.2 Principle of Operation	10
3.2 Reference to previous generation(s) or variants.....	10
3.3 Accessories Details	11
3.4 Combination with other Medical Devices	11
4. Risks and warnings	11
4.1 Residual risks and undesirable side effects	11

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	4 of 31

4.2	Warnings.....	11
4.3	Precautions.....	12
4.4	Other relevant aspects of safety	12
5	Summary of clinical evaluation and post-market clinical follow-up (PMCF)	12
5.1	Summary of clinical data related to similar device, if applicable	12
5.3	Summary Of Clinical Data from Conducted Investigations of the device before The CE-Marking 13	
5.4	Summary of clinical data from other sources	13
5.5	An overall summary of the clinical performance and safety.....	16
5.6.	Ongoing or planned post-market clinical follow-up.....	17
6.	Possible diagnostic or therapeutic alternative.....	18
7.	Suggested profile and training for users	18
8.	Reference to any harmonized standards and CS applied	19
8.1	Applicable Harmonized Standards	19
8.2	Other Applicable Standards.....	20
8.3	List of Guidelines	22
	Summary Of Safety and Clinical Performance	24
	Intended for Patients	24
	Introduction.....	24
1.	Device identification and general information	24
2.	Intended use of the device.....	25
2.1	Intended Purpose	25
2.2	Indications & Target Populations	25
2.3	Contraindications	25
3.	Device Description.....	26
3.1	Device description	26
3.2	Materials that come in contact with patient.....	26

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	5 of 31

3.3	Information about medicinal substances in the device, if any.....	26
3.4	Description of how the device is achieving its intended mode of action	27
3.5	Description of accessories, if any	27
4.	Risks and warnings	27
4.1	How potential risks have been controlled or managed	27
4.2	Residual Risks.....	27
4.3	Adverse events	28
4.4	Warnings.....	28
4.5	Precautions	28
4.6	Summary of any field safety corrective action, (FSCA including FSN) if applicable	29
5.	Summary of clinical evaluation and post-market clinical follow-up	29
5.1	Clinical background of the device.....	29
5.2	The clinical evidence for the CE-marking	29
5.3	Safety	30
9.	Possible diagnostic or therapeutic alternative.....	30
10.	Suggested profile and training for users	31
11.	Revision history	31

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	6 of 31

Introduction

The Regulation (EU) 2017/745 on medical devices requires that the manufacturer shall draw up a summary of safety and clinical performance (SSCP) for implantable devices and class III devices, other than custom-made or investigational devices. The SSCP shall be validated by a notified body (NB) and made available to the public via the European database on medical devices (Eudamed).

The SSCP is intended to provide public access to an updated summary of clinical data and other information about the safety and clinical performance of the medical device. The SSCP will be an important source of information for intended users – both healthcare professionals and relevant for patients. It is one of several means intended to fulfill the objectives of the Medical Device Regulation (MDR) to enhance transparency and provide adequate access to information.

The SSCP is not intended to:

- Give general advice on the diagnosis or treatment of particular medical conditions, nor
- Replace the instructions for use (IFU) as the main document that will be provided to ensure the safe use of a particular device, nor
- Replace the mandatory information on implant cards or in any other mandatory documents.

The main purpose of this document is to guide the presentation, content and validation of the SSCP. The word “shall” is used when there is a corresponding “shall” in the MDR, otherwise “should” or “recommended” etc. is used to indicate the interpretation of the MDR.

The following information is intended for users/healthcare professionals.

1. Device identification and general information

1.1 Device Trade Name(s)

Product Name	Urethral Stents-With/without Hydrophilic coated
Brand Name:	Devon
Variant Name:	-----

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	7 of 31

1.2 Manufacturer's Name & Address

Legal Manufacturer Name:	DEVON INNOVATIONS PRIVATE LIMITED
Registered Office & Manufacturing Unit-I Address:	No. 27A, Near State Bank of India, Electronic City Phase I, Hosur Main Road, Bangalore-560 100, India. Phone no: 080-28522354/28522367/28522368
Manufacturing Unit-II Address:	Gupta complex, 1st floor, Khasra No: 519/370, Near EWS flats, sector-1, village Kamli Parwanoo 173220 Himachal Pradesh, India. Phone no: 01792232492
Email:	srinivas@devoncath.com , nagendrakumar@devoncath.com
Website:	www.devoncath.com

1.3 Manufacturer's Single Registration Number (SRN)

Single Registration Number (SRN) for Unit-I Manufacturing site:	IN-MF-000010584
Single Registration Number (SRN) for Unit-II Manufacturing site:	IN-MF-000045808

1.4 Basic UDI-DI

8903410USW6

1.5 Medical Device Nomenclature

EMDN Code:	U020399
MDN Code:	MDN 1104-2
MDS and MDT code:	MDS 1005, MDT 2001, MDT 2002, MDT 2008, MDT 2011

1.6 Class of device

Class IIb, Rule 05, in accordance with Annex VIII of EU Medical Device Regulation 2017/745

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	8 of 31

1.7 Year of first certificate (CE) of the subject device

2012

1.8 Authorized Representative

Name:	Amstermed BV
Address:	Saturnusstraat 46-62, Unit 032,2132 HB Hoofddorp The Netherlands.
Phone:	+31 23 56 56 337
Email:	info@amstermed.nl
Website:	https://www.amstermed.nl
SRN:	NL-AR-000001971

1.9 NB Details

Name:	DNV Product Assurance AS
Address:	Veritasveien 1, 1363 Høvik, Norway
Website:	www.dnv.com
Notified Body No.:	2460

1.10 Conformity Assessment Procedure

Conformity assessment procedure followed is Annexure IX.

1.11 Link to SSCP in website

The link for the Summary of Safety and Clinical Performance (SSCP) is provided below:

2. Intended use of the device

2.1 Intended Purpose

Used for stenting the Urethra after Hypospadias or Epispadias repair and to provide postoperative drainage of the bladder. In case of Hydrophilic coated, it is to improve the ease of insertion.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	9 of 31

2.2 Indications & Target Populations

Indication:

- Hypospadias or Epispadias repair
- Provide postoperative drainage of the bladder

Target patient population: Adult

2.3 Contraindications

- Contraindicated surgical candidate
- Unexplained hematuria
- Unrepaired ureteral avulsion

3. Device Description

3.1 Description of the Device

A stent is a hollow tube that maintains patency until healing can take place or an obstruction is relieved.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	10 of 31

Product Image:

Product Name	Packing	Image
Urethral Stents-With/without Hydrophilic coated	With Packing	
	Without Packing	

3.1.1 Device Models

Not Applicable. The Urethral Stents- With/Without Hydrophilic coated does not have any variants.

3.1.2 Principle of Operation

A Urethral Stents-With/without Hydrophilic coated is a small flexible tube inserted into the urethra to keep it open and allow for proper urine flow, to treat urinary tract obstruction or blockage.

3.2 Reference to previous generation(s) or variants

Legacy Device Name:	Stents (Urology)- Urethral Stents- With/Without Hydrophilic coated
Brand/Proprietary Name:	Devon
93/42/EEC (MDD) Cert. No.:	246182-2017-CE-IND-NA-PS, Rev.2.0
Notified Body Details:	DNV Product Assurance AS
Is any significant difference between Legacy Device & Device Under Evaluation?	There is no significant difference between legacy device and subject device with respect to raw materials used in production, device description, intended purpose, medical indications, target user, target patient population, side-effects, and contraindications.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	11 of 31

3.3 Accessories Details

Not Applicable. The Urethral Stents- With/without Hydrophilic coated does not have any accessories supplied by the manufacturer.

3.4 Combination with other Medical Devices

Urethral Stents- With/without Hydrophilic coated is used along with Guidewire. The stent can be used with any compatible guidewire from other manufacturers.

4. Risks and warnings

4.1 Residual risks and undesirable side effects

a. Residual Risks

- Urinary Tract Infection
- Stent migration
- Contamination or Deterioration of product
- Toxic to environment
- Stent fragmentation
- Haematuria
- Tissue damage
- Delay in procedure, Inconvenience to the user

b. Adverse Events

- Migration
- Sepsis
- Encrustation

4.2 Warnings

- Urethral Stents are intended for single-use only.
- Duration of Use:
 - Periodic evaluation is advised. The Stent must not remain indwelling more than three months. These stents are not indented as permanent indwelling devices

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	12 of 31

- Do not use device if there is any indication that the sterility of the device has been compromised.
- Adverse effects
- Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance. Follow up procedures.
- Reuse:
- Reusing single-use stents can lead to urinary tract infections in patients.

4.3 Precautions

Carefully read all instructions for use and product labeling. The device shall only be applied for its intended use, and in accordance with these instructions. Observe all cautions and warnings throughout these instructions. Failure to do so may result in complications.

All Health care professionals is responsible for using the appropriate technique and deciding on the indication for use of this device based on own experience, training and medical judgment. The doctor must be trained in the proper use of the device.

4.4 Other relevant aspects of safety

There were no identified and/or received reportable events that led to death, a serious deterioration in the state of health of the patient, user, or other person for Urethral Stents- With/without Hydrophilic coated. Hence FSCA or FSN is not applicable.

5 Summary of clinical evaluation and post-market clinical follow-up (PMCF)

5.1 Summary of clinical data related to similar device, if applicable

The Urethral Stents- With/without Hydrophilic coated belongs to the “Urology” group. In the present market there are many similar devices and/or benchmark devices available with same intended purpose and are having the same generally acknowledged state-of-the-art.

These similar devices fall under “Well-Established Technology”. Data from similar devices is considered for the conformation of conformity to the Urethral Stents- With/without Hydrophilic coated relevant general safety and performance requirements. The similar device data is used to demonstrate ubiquity of design, lack of novelty, known safety and performance profile of a generic group of devices, etc.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	13 of 31

The below mentioned similar devices data is used to evaluate the Urethral Stents- With/without Hydrophilic coated relevant general safety and performance requirements as part of literature review. These similar devices contain the same Raw materials and same intended purpose, but due to insufficient information availability the full assessment of equivalence is not possible. Therefore, these similar devices are used for the same clinical intended purposes as the Urethral Stents- With/without Hydrophilic coated and are considered to be similar but non-equivalent devices.

S.NO	Product Name	Variant Name	Similar Device	Manufacturer Name
1.	Urethral Stents- With/without Hydrophilic coated	-----	Zaontz Urethral Stent	Cook® Medical

The similar device's SSCP would be available in EUDAMED.

5.3 Summary Of Clinical Data from Conducted Investigations of the device before The CE-Marking

Not Applicable

5.4 Summary of clinical data from other sources

The below mentioned are literatures selected for detailed review for:

- Evaluation of state of the art
- Evaluation of clinical data from similar devices

#	ID#	Source Link	Literature Title
1.	L1	https://pubmed.ncbi.nlm.nih.gov/14688404/	Ureteral Stent Placement without Postprocedural Nephrostomy Tube: Experience in 41 Patients
2.	L3	https://pubmed.ncbi.nlm.nih.gov/32354034/	Novel Treatment Strategy for Management of Traumatic Bulbar Urethral Rupture Using Temporary Urethral Stent after Primary Realignment; Retrospective Comparison between Thermo-Expandable Urethral Stent and Self-Expandable Polymer-Coated Urethral Stent

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	14 of 31

#	ID#	Source Link	Literature Title
3.	L5	https://pubmed.ncbi.nlm.nih.gov/16646656/	Comparison of Nitinol Urethral Stent Infections with Indwelling Catheter-Associated Urinary-Tract Infections
4.	L7	https://pubmed.ncbi.nlm.nih.gov/12629371/	Experience with urethral stent explantation
5.	L9	https://pubmed.ncbi.nlm.nih.gov/32293992/	Clinical Effectiveness of Single Pigtail Suture Stent on Patient Comfort: A Double-Blind Prospective Randomized Trial
6.	L11	https://pubmed.ncbi.nlm.nih.gov/31549458/	Ureteric stents: Overview of current clinical applications and economic implications
7.	L12	https://pubmed.ncbi.nlm.nih.gov/29234836/	Trans-Urethral Ureteral Stent Replacement Technique (TRUST): 10-Year Experience in 1168 Patients
8.	L13	https://www.academia.edu/download/80791286/3_mayank-.pdf	A prospective study comparing stenting versus non-stenting after uncomplicated ureteroscopy for ureteric stones
9.	L14	https://pubmed.ncbi.nlm.nih.gov/34245667/	A randomized trial investigating clinical outcomes and stent-related symptoms after placement of a complete intra-ureteric stent on a string versus conventional stent placement
10.	L15	https://pubmed.ncbi.nlm.nih.gov/24419714/	Memokath™ Urethral Stents Induce Incontinence in Patients with Urethral Balloon Catheters
11.	L16	https://pubmed.ncbi.nlm.nih.gov/27256191/	The long-term results of temporary urethral stent placement for the treatment of recurrent bulbar urethral stricture disease
12.	L17	https://pubmed.ncbi.nlm.nih.gov/29751026/	Urethroplasty after Urethral Urolume Stent: an International Multicenter Experience
13.	L18	https://pubmed.ncbi.nlm.nih.gov/17937962/	Preservation of Lower Urinary Tract Function in Posterior Urethral Stenosis: Selection of Appropriate Patients for Urethral Stents
14.	L19	https://pubmed.ncbi.nlm.nih.gov/25419377/	Management of recurrent bulbar urethral stricture-a 54 patients' study with Allium bulbar urethral stent (BUS)

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	15 of 31

Examples of clinical data registries are:

1. Clinicaltrials.gov: ClinicalTrials.gov is a database of privately and publicly funded clinical studies conducted around the world. ClinicalTrials.gov is a resource provided by the U.S. National Library of Medicine.
2. Clinicaltrialsregister.eu: The European Union Clinical Trials Register allows you to search for protocol and results information on interventional clinical trials that are conducted in the European Union (EU) and the European Economic Area (EEA); clinical trials conducted outside the EU / EEA that are linked to European pediatric-medicine development.
3. Ctri.gov: The Clinical Trials Registry- India (CTRI), hosted at the ICMR's National Institute of Medical Statistics (<http://icmr-nims.nic.in>), is a free and online public record system for registration of clinical trials being conducted in India. Any researcher who plans to conduct a trial involving human participants, of any intervention such as drugs, surgical procedures, preventive measures, lifestyle modifications, devices, educational or behavioral treatment, rehabilitation strategies as well as trials are expected to register the trial in the CTRI before enrollment of the first participant. In the CTRI, details of Indian investigators, trial sites, Indian target sample size and date of enrollment are captured.

Studies referred from device registries are given below:

Sl.no	Registry	Study No	Country	Study Title	Study Phase	Type of study	Number of subjects
1.	https://clinicaltrials.gov/study/NCT02323243?term=urethra&limit=50&cond=stent&rank=1	NCT02323243	France	Clinical Interest of the Temporary Urethral Stent Allium " Bulbar Urethral Stent " in the Treatment of Detrusor Sphincter Dyssynergia	Phase1 Phase 2	Interventional	9
2.	https://clinicaltrials.gov/study/NCT00252941?term=urethra&limit=50&cond=stent&rank=3	NCT00252941	United States	Prophylactic Urethral Stenting With Memokath After Prostate Implantation for Prostate Adenocarcinoma	Phase1 Phase 2	Interventional	20

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	16 of 31

--	--	--	--	--	--	--	--

The advantage of referring the data from a third-party registry is usually centered on patients or cases, which provides a sufficient level of traceability and detail for the study of our product. Using registries can be less expensive than initiating new and private device registries. The data in these registries might be relevant but may not contain the exact attributes needed to document/demonstrate clinical performance.

Hence, we cannot completely rely on the above registries to collect enough data for our devices. We have referred these registries only to gather information on similar devices and the study methodology.

5.5 An overall summary of the clinical performance and safety

Residual Risks	Medical Benefits
1. Urinary Tract Infection 2. Stent migration 3. Contamination or Deterioration of product 4. Toxic to environment 5. Stent fragmentation 6. Haematuria 7. Tissue damage 8. Delay in procedure, Inconvenience to the user	Subject Device 1. Relief of Urinary Obstruction 2. Continuous urine drainage 3. Prevention of Ureteral Stricture Formation Similar Device 1. Stent removal is technically feasible. 2. Minimal local tissue reaction. 3. Improved QoL and less requirement of analgesic 4. lower ureteral irritation and decreases VUR (Vesicoureteral reflux) 5. Improved pain and general health 6. Minimises postoperative ureteral obstruction 7. Promote healing 8. Reduce the incidence of ureteral stricture 9. Decreases distal migration of stents Reduces SRSs (Stent related symptoms)

All the residual risks were reviewed and analysed and also identified the medical benefits of the intended use outweigh the overall residual risk. Urethral Stents-With/without Hydrophilic coated complies with all the Safety and Performance requirements with respect to the intended use of the device from the Clinical Evaluation study. The clinical evaluation is complete and conforms to the general safety and performance requirements. The Clinical evidence is demonstrated with the relevant General Safety & Performance

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	17 of 31

Requirements as per Annex I of MDR. Risk Mitigation has been established as per the guidelines of ISO 14971.

5.6. Ongoing or planned post-market clinical follow-up

Overall, 35 subjects have participated in this PMCF study of the product Urethral Stent With/without Hydrophilic coated. As per the PMCF study, the Urethral Stent With/without Hydrophilic coated have met the primary and secondary objectives.

Overall study results:

Parameter	Study Results												
Subjects	A total of 35 Subjects were enrolled for the PMCF study. As per the study, 30 (86%) subjects were male and 5 (14%) subjects were female.												
	Subjects Group summary –35 (100%) subjects were adults.												
	Age Group:												
	<table><tr><th>Age Group</th><th>No: of subjects</th><th>Percentage (%)</th></tr><tr><td>31-40 Years</td><td>12</td><td>34%</td></tr><tr><td>41-50 Years</td><td>19</td><td>54%</td></tr><tr><td>51-60 Years</td><td>4</td><td>11%</td></tr></table>	Age Group	No: of subjects	Percentage (%)	31-40 Years	12	34%	41-50 Years	19	54%	51-60 Years	4	11%
	Age Group	No: of subjects	Percentage (%)										
31-40 Years	12	34%											
41-50 Years	19	54%											
51-60 Years	4	11%											
Target Users	Urologist												
Study Site	Kapoor’s Kidney Centre Pvt Ltd, Chandigarh, National Kidney hospital Jalandhar and Yasoda hospital Pune												
Clinical Indication	Provide postoperative drainage of the bladder and Hypospadias or Epispadias repair												
Clinical Safety	All the subjects considered for this study benefitted out of using Urethral Stent with/without hydrophilic coating. This clinically proves the safety of using our product on subjects with better efficacy.												
Clinical Performance	The overall rating was 4 which is “Good” as per the definition, hence the product performance was clinically proven to be “Good”.												
Adverse events and Risks	No users reported any of the adverse events during the study.												
	No reports on the risks imposed by our Urethral Stent With/without Hydrophilic coated.												

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	18 of 31

Follow Up Summary	The stent was removed during a follow up procedure. The subject had no issues with stent removal. The Urethral Stent with/without hydrophilic coating effectively relieved the obstructive symptoms, and the patient experienced significant symptom improvement after the procedure. No complications were observed during follow-up.
Product Experience	The overall rating is 8, it is “Good” as per the definition. Hence this proves that our users were satisfied with the product and its purpose.

The Urethral Stent With/without Hydrophilic coated from Devon Innovations Private Limited has reached all the safety and performance requirements with respect to the intended use of the device from the Post Market Follow up Clinical study. The Performance and Safety of the device as we claimed have been established in the Technical File and Instruction for Use (IFU). There were no new risks identified from the PMCF study for the product hence there is no addition to the residual risks which we have already identified in the Risk Management Report and that is been mitigated and are acceptable when weighed against the benefits to the patient.

6. Possible diagnostic or therapeutic alternative

- Urethroplasty
- Direct Visual Internal Urethrotomy (DVIU)
- Urethral Dilation

7. Suggested profile and training for users

Target Users: Urologist

Urologists diagnose and treat diseases of the urinary tract in both men and women. They also diagnose and treat anything involving the reproductive tract in men.

The target users are aware of basic operations of Urethral Stent With/without Hydrophilic coated. There is no special user training is required. However, the device related directions for use information are provided in the Instruction for Use.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	19 of 31

8. Reference to any harmonized standards and CS applied

8.1 Applicable Harmonized Standards

#	Standard ID	Current Issue	Title
Quality Management System Requirements			
1.	EN ISO 13485	2016/AC:2018/A11:2021	Medical devices - Quality management systems - Requirements for regulatory purposes (ISO 13485:2016)
Risk Management Requirements			
2.	EN ISO 14971	2019/A11:2021	Medical devices - Application of risk management to medical devices (ISO 14791:2019)
Biological Risk Evaluation Requirements			
3.	EN ISO 10993-10	2023	Biological evaluation of medical devices - Part 10: Tests for skin sensitization (ISO 10993-10:2021)
4.	EN ISO 10993-23	2021	Biological evaluation of medical devices - Part 23: Tests for irritation (ISO 10993-23:2021)
Labels & Symbols Requirements			
5.	EN ISO 15223-1	2021	Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements
Packaging Requirements			
6.	EN ISO 11607-1	2020/A1:2023	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems - (ISO 11607-1:2019)
7.	EN ISO 11607-2	2020/A1:2023	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes (ISO 11607-2:2019)
EO Sterilization Requirements			
8.	EN ISO 11135	2014/A1:2019	Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices (ISO11135:2014)
Sterility Test Requirements			

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	20 of 31

#	Standard ID	Current Issue	Title
9.	EN ISO 11737-1	2018/A1:2021	Sterilization of health care products - Microbiological methods - Part 1: Determination of a population of microorganisms on products - (ISO 11737-1:2018)
10.	EN ISO 11737-2	2020	Sterilization of health care products - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process (ISO 11737-2:2019)

8.2 Other Applicable Standards

#	Standard ID	Current Issue	Title
Risk Management Requirements			
1.	ISO/TR 24971	2020	Medical devices – Guidance on the application of ISO 14971 (ISO/TR 24971:2020)
Usability Requirements			
2.	EN 62366-1	2015/A1:2020	Medical devices - Part 1: Application of usability engineering to medical devices
3.	IEC 62366-1	2015/A1:2020	Medical devices — Part 1: Application of usability engineering to medical devices Amendment 1
Biological Risk Evaluation Requirements			
4.	EN ISO 10993-1	2020	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process (ISO 10993-1:2018, including corrected version 2018-10)
5.	EN ISO 10993-3	2014	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity (ISO 10993-3:2014)
6.	EN ISO 10993-5	2009/A11:2025	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
7.	EN ISO 10993-6	2016	Biological evaluation of medical devices - Part 6: Tests for local effects after implantation (ISO 10993-6:2016)

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	21 of 31

#	Standard ID	Current Issue	Title
8.	EN ISO 10993-7	2008/A1:2022	Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals - Amendment 1: Applicability of allowable limits for neonates and infants (ISO 10993-7:2008/Amd 1:2019)
9.	EN ISO 10993-11	2018	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity (ISO 10993-11:2017)
Instructions For Use Requirements			
10.	EN ISO 20417	2021	Medical devices - Information to be supplied by the manufacturer (ISO 20417:2021, Corrected version 2021-12)
Medical Device "Sterile" Requirements			
11.	EN 556-1	2024	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
Transport Requirements			
12.	ISTA 2A	2011	Partial Stimulation Performance Test Procedure Packaged Products weighing 150 lbs (68 kg) or less
Cleanroom Requirements			
13.	ISO 14644-1	2015	Clean Rooms and Associated Controlled Environments – Part 1: Classification of Air Cleanliness by Particle Concentration
14.	ISO 14644-2	2015	Clean Rooms and Associated Controlled Environments – Part 2: Monitoring to Provide Evidence of Clean Room Performance Related to Air Cleanliness by Particle Concentration
Post Market Surveillance Requirements			
15.	ISO/ TR 20416	2020	Medical devices – Post market surveillance for manufacturers
Stability Requirements			
16.	ASTM F 1980-21	2021	Standard Guide for Accelerated Aging of Sterile Barrier Systems and Medical Devices
Labels & Symbols Requirements			

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	22 of 31

#	Standard ID	Current Issue	Title
17.	ISO 15223-1	2021/Amd 1:2025	Medical devices — Symbols to be used with information to be supplied by the manufacturer — Part 1: General requirements Amendment 1: Addition of defined term for authorized representative and modified EC REP symbol to not be country or region specific

8.3 List of Guidelines

#	Guideline	Current Issue	Title
1.	MEDDEV 2.7.1 Rev. 4	June 2016	Clinical Evaluation: A guide for manufacturers and notified bodies under directives 93/42/EEC and 90/385/EEC
2.	MEDDEV 2.5/5 Rev. 3	February 1998	Translation Procedure - Guidelines relating to the application of: The council directive 90/385/EEC on active implantable medical devices The council directive 93/42/EEC on medical devices
3.	MEDDEV 2.12-1 Rev. 8	January 2013	Guidance document - Market surveillance - Guidelines on a Medical Devices Vigilance System
4.	NB-MED 2.12/Rec. 1	February 2000	Post-Marketing Surveillance (PMS) post market/production
5.	MDCG 2021-24	Oct 2021	Guidance on classification of medical devices
6.	MDCG 2020-6	2020	Clinical evidence needed for medical devices previously CE marked under Directives 93/42/EEC.
7.	MDCG 2018-1 Rev.4	April 2021	Guidance on BASIC UDI-DI and changes to UDI-DI
8.	MDCG 2020-7	April 2020	Post-market clinical follow-up (PMCF) Plan Template - A guide for manufacturers and notified bodies
9.	MDCG 2020-8	April 2020	Post-market clinical follow-up (PMCF) Evaluation Report Template - A guide for manufacturers and notified bodies

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	23 of 31

#	Guideline	Current Issue	Title
10.	MDCG 2022-21	December 2022	Guidance on Periodic Safety Update Report (PSUR) according to Regulation (EU) 2017/745 (MDR)
11.	MDCG 2019-9 Rev.1	March 2022	Summary of safety and clinical performance A guide for manufacturers and notified bodies
12.	MDCG 2020-3 Rev.1 6	September 2023	Guidance on significant changes regarding the transitional provision under Article 120 of the MDR - May 2023
13.	MDCG 2024-2	Feb 2024	MDCG 2024-2 Procedures for the updates of the European Medical Device Nomenclature

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	24 of 31

Summary Of Safety and Clinical Performance

Intended for Patients

Introduction

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of its safety and clinical performance is prepared for users/healthcare professionals.

The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your healthcare professional in case you have questions about your medical condition or about the use of the device in your situation. This SSCP is not intended to replace an Implant card or the Instructions for Use to provide information on the safe use of the device.

The following information is intended for patient.

READABILITY

The SSCP has one part for intended users/healthcare professionals, and a second part for patients. It provides clear information at an appropriate depth to reflect the healthcare professionals' and the patients' different levels of knowledge. The readability of the part of the SSCP intended for patients is assessed by a test given to lay persons.

The readability was assessed by providing the document to laypersons with different literacy background to read. All the laypersons who have read the document were able to understand the terms and details of the device properly. Also, they agreed that the document language was legible and understandable. The manufacturer used this method to confirm that the SSCP was written in a way that is clear to the patient.

1. Device identification and general information

Product Name	Urethral Stents-With/without Hydrophilic coated
Brand Name:	Devon
Models:	-----

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	25 of 31

Manufacturer's Name & Address	DEVON INNOVATIONS PRIVATE LIMITED Registered Office & Manufacturing Unit-I Address: No. 27A, Near State Bank of India, Electronic City Phase I, Hosur Main Road, Bangalore-560 100, India. Phone no: 080-28522354/28522367/28522368 Manufacturing Unit-II Address: Gupta complex, 1st floor, Khasra No: 519/370, Near EWS flats, sector-1, village Kamli Parwanoo 173220 Himachal Pradesh, India. Phone no: 01792232492 Email: srinivas@devoncath.com nagendrakumar@devoncath.com Website: www.devoncath.com
Basic UDI-DI	8903410USW6
Year when the device was first CE-marked	2012

2. Intended use of the device

2.1 Intended Purpose

Used for stenting the Urethra after Hypospadias or Epispadias repair and to provide postoperative drainage of the bladder. In case of Hydrophilic coated, it is to improve the ease of insertion.

2.2 Indications & Target Populations

Indication:

- Hypospadias or Epispadias repair
- Provide postoperative drainage of the bladder

Target patient population: Adult

2.3 Contraindications

- Contraindicated surgical candidate
- Unexplained hematuria

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	26 of 31

- Unrepaired ureteral avulsion

3. Device Description

3.1 Device description

A stent is a hollow tube that maintains patency until healing can take place or an obstruction is relieved.

Product Image:

Product Name	Packing	Image
Urethral Stents- With/without Hydrophilic coated	With Packing	
	Without Packing	

3.2 Materials that come in contact with patient

Material that comes in contact with patient is stent.

3.3 Information about medicinal substances in the device, if any

Urethral Stents-With/without Hydrophilic coated does not incorporate medicinal substances. Hence this declaration is not applicable

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	27 of 31

3.4 Description of how the device is achieving its intended mode of action

Urethral stent is placed in the urethra to prevent or relieve obstruction. These stents are deployed to keep the urethral passage open and allow normal urinary flow, particularly after surgeries or in cases of urethral strictures.

3.5 Description of accessories, if any

Not Applicable. The Urethral Stents- With/without Hydrophilic coated does not have any accessories supplied by the manufacturer

4. Risks and warnings

Contact your healthcare professional if you believe that you are experiencing side-effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

4.1 How potential risks have been controlled or managed

We have established and maintaining a strong quality management system and continuous process monitoring controls as per EN ISO 13485:2016/A11:2021 requirements. As per the risk management plan, all the necessary and possible risk control measures are implemented to reduce the risk to practicable acceptable level. We have implemented all necessary the risk control measures for each hazard for reducing the risks to an acceptable level.

The Urethral Stents- With/without Hydrophilic coated risk management is completed by risk estimation, risk analysis, risk evaluation, risk control, overall residual risk evaluation, production, and postproduction information. The Urethral Stents- With/without Hydrophilic coated is reached all the safety and performance claims when used as per the defined indented purpose.

The Urethral Stents- With/without Hydrophilic coated medical benefits and residual risks are compared based on medical condition, literature data certainty, similar device data (benefits & residual risks) from the market, acknowledged state of the art.

4.2 Residual Risks

- Urinary Tract Infection

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	28 of 31

- Stent migration
- Contamination or Deterioration of product
- Toxic to environment
- Stent fragmentation
- Haematuria
- Tissue damage
- Delay in procedure, Inconvenience to the user

4.3 Adverse events

- Migration
- Sepsis
- Encrustation

4.4 Warnings

- Urethral Stents are intended for single-use only.
- Duration of Use:
 - Periodic evaluation is advised. The Stent must not remain indwelling more than three months. These stents are not intended as permanent indwelling devices
 - Do not use device if there is any indication that the sterility of the device has been compromised.
- Adverse effects
 - Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance. Follow up procedures.
- Reuse:
 - Reusing single-use stents can lead to urinary tract infections in patients.

4.5 Precautions

Carefully read all instructions for use and product labeling. The device shall only be applied for its intended use, and in accordance with these instructions. Observe all cautions and warnings throughout these instructions. Failure to do so may result in complications.

All Health care professionals is responsible for using the appropriate technique and deciding on the indication for use of this device based on own experience, training and medical judgment. The doctor must be trained in the proper use of the device.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	29 of 31

4.6 Summary of any field safety corrective action, (FSCA including FSN) if applicable

There were no identified and/or received reportable events that led to death, a serious deterioration in the state of health of the patient or user, for Urethral Stents-With/without Hydrophilic coated. Hence FSCA or FSN is not applicable.

5. Summary of clinical evaluation and post-market clinical follow-up

5.1 Clinical background of the device

The Urethral Stents- With/Without Hydrophilic coated is a designed and developed as per the latest and/or current technical, international and regulatory. The Urethral Stents- With/Without Hydrophilic coated performance complies all necessary requirements required for its intended purpose. The Urethral Stents- With/Without Hydrophilic coated related all known foreseeable hazards are identified and associated risk are reduced as far as possible by implementing all necessary risk control measures as per the requirements of EN ISO 14971:2019/A11:2021. The Urethral Stents- With/Without Hydrophilic coated is manufactured as per the defined standard operating procedures in controlled environments by training personnel and complies all necessary requirements of EN ISO 13485:2016/A11:2021.

5.2 The clinical evidence for the CE-marking

Urethral Stents- With/Without Hydrophilic coated is having a CE certified medical device under EU MDD 93/42/EEC. Hence, the Urethral Stents- With/Without Hydrophilic coated comply the definition as a legacy device as per the MDCG 2020-6:2020 – Regulation (EU) 2017/745: Clinical evidence needed for medical devices previously CE marked under Directives 93/42/EEC or 90/385/EEC.

Legacy Device Name:	Stents (Urology)
Brand/Proprietary Name:	Devon
93/42/EEC (MDD) Cert. No.:	246182-2017-CE-IND-NA-PS, Rev.2.0
Notified Body Details:	DNV Product Assurance AS
Is any significant difference between Legacy Device & Device Under Evaluation?	There is no significant difference between legacy device and subject device with respect to raw materials used in production, device description, intended purpose, medical indications, target user, target patient population, side-effects, and contraindications.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	30 of 31

The Urethral Stents- With/Without Hydrophilic coated belongs to the “Urology” group. In the present market there are many similar devices and/or benchmark devices available with same intended purpose and are having the same generally acknowledged state-of-the-art.

S.NO	Product Name	Variant Name	Similar Device	Manufacturer Name
1.	Urethral Stents- With/without Hydrophilic coated	-----	Zaontz Urethral Stent	Cook® Medical

5.3 Safety

We have reviewed and analysed all the residual risks and also identified the medical benefits of the intended use outweigh the overall residual risk. All the residual risks are acceptable by providing appropriate information to the end user’s awareness in the form of “Label’ and ‘Instructions for Use”.

Residual Risks	Medical Benefits
1. Urinary Tract Infection 2. Stent migration 3. Contamination or Deterioration of product 4. Toxic to environment 5. Stent fragmentation 6. Haematuria 7. Tissue damage 8. Delay in procedure, Inconvenience to the user	Subject Device 1. Relief of Urinary Obstruction 2. Continuous urine drainage 3. Prevention of Ureteral Stricture Formation Similar Device 1. Stent removal is technically feasible. 2. Minimal local tissue reaction. 3. Improved QoL and less requirement of analgesic 4. lower ureteral irritation and decreases VUR (Vesicoureteral reflux) 5. Improved pain and general health 6. Minimises postoperative ureteral obstruction 7. Promote healing 8. Reduce the incidence of ureteral stricture 9. Decreases distal migration of stents 10. Reduces SRSs (Stent related symptoms)

9. Possible diagnostic or therapeutic alternative

When considering alternative treatments, it is recommended to contact your healthcare professional who can take into account your individual situation.

	Devon Innovations Private Limited	Document No.:	TD/DIP/SSCP/02
		Revision No.:	00
	Urethral Stents-With/without Hydrophilic coated	Effective Date:	23.06.2025
	Summary of Safety and Clinical Performance	Page No.:	31 of 31

10. Suggested profile and training for users

Target Users: Urologist

Urologists diagnose and treat diseases of the urinary tract in both men and women. They also diagnose and treat anything involving the reproductive tract in men.

The target users are aware of basic operations of Urethral Stents-With/without Hydrophilic coated. There is no special user training is required. However, the device related directions for use information are provided in the Instruction for Use.

11. Revision history

SSCP Rev. No.	Date Issued	Change description	Rev. Validated by the NB
00	23.06.2025	Initial Release	☐ Yes Validation language: English ☒ No